



TRIM3 Polyclonal Antibody

Catalog No	YP-Ab-04255
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	TRIM3
Immunogen	The antiserum was produced against synthesized peptide derived from human TRIM3. AA range:1-50
Specificity	TRIM3 Polyclonal Antibody detects endogenous levels of TRIM3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNF22;RING finger protein 97;RING finger protein 22;HAC1;FLJ16135;TRIM3; Tripartite motif-containing protein 3;BERP;RNF97;Brain-expressed RING finger protein
Observed Band	81kD
Calculated Molecular Weight	81kD
Cell Pathway	Cytoplasm . Early endosome . Golgi apparatus, trans-Golgi network . Cell projection, dendrite .
Tissue Specificity	Expressed in brain, heart, uterus and testis.
Function	domain:The interaction with myosin V is dependent upon its NHL repeats, which form a beta-propeller (NHL) domain containing six blades.,similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger., similarity:Contains 1 filamin repeat.,similarity:Contains 1 RING-type zinc finger., similarity:Contains 6 NHL repeats.,subunit:Associates with myosin V and alpha-actinin-4 (ACTN4).,tissue specificity:Expressed in brain, heart, uterus and testis.,



Background

The protein encoded by this gene is a member of the tripartite motif (TRIM) family, also called the 'RING-B-box-coiled-coil' (RBCC) subgroup of RING finger proteins. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to cytoplasmic filaments. It is similar to a rat protein which is a specific partner for the tail domain of myosin V, a class of myosins which are involved in the targeted transport of organelles. The rat protein can also interact with alpha-actinin-4. Thus it is suggested that this human protein may play a role in myosin V-mediated cargo transport. Alternatively spliced transcript variants encoding the same isoform have been identified. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images




Western Blot analysis of COLO cells using TRIM3 Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000




Immunofluorescence analysis of NIH/3T3 cells, using TRIM3 Antibody. The picture on the right is blocked with the synthesized peptide.

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TRIM3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from COLO cells, using TRIM3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from Jurkat cells using TRIM3 antibody.