



FA5 (heavy chain, Cleaved-Ala29) rabbit pAb

Catalog No	YP-Ab-04373
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	F5
Immunogen	Synthesized peptide derived from human FA5 (heavy chain, Cleaved-Ala29)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat FA5 (heavy chain, Cleaved-Ala29, protein was cleaved amino acid sequence between 28-29)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA5;F5;FA5 (heavy chain;Cleaved-Ala29);Coagulation factor V (Activated protein C cofactor;Proaccelerin;labile factor) [Cleaved into: Coagulation factor V heavy chain;Coagulation factor V light chain]
Observed Band	245kD,80kD
Calculated Molecular Weight	245kD,80kD
Cell Pathway	Secreted .
Tissue Specificity	Plasma.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
