



NOTCH3 (Cleaved-Val1662) rabbit pAb

Catalog No	YP-Ab-04397
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	NOTCH3
Immunogen	Synthesized peptide derived from human NOTCH3 (Cleaved-Val1662)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat NOTCH3 (Cleaved-Val1662, protein was cleaved amino acid sequence between 1661-1662)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NOTCH3;NOTCH3 (Cleaved-Val1662);Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation;Notch 3 intracellular domain]
Observed Band	75 260kD
Calculated Molecular Weight	75 260kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein.; [Notch 3 intracellular domain]: Nucleus. Following proteolytical processing NICD is translocated to the nucleus.
Tissue Specificity	Ubiquitously expressed in fetal and adult tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images