



FLAG-tag Polyclonal Antibody

Catalog No	YP-Ab-04401
Isotype	IgG
Reactivity	Species independent
Applications	WB;ELISA
Gene Name	Flag Tag
Immunogen	DDDDK synthetic peptide conjugated to KLH.
Specificity	FLAG-tag Polyclonal Antibody detects FLAG-tagged recombinant proteins or FLAG-tagged proteins overexpressed in cells.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.05% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/1000 - 1/3000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Flag Tag;Flag-Tag DDDDK TAG, DDDDK-TAG;DDDDK epitope tag; DDDDK epitope tag; DYKDDDDK epitope tag
Observed Band	55KDa
Calculated Molecular Weight	55KD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
