




BSA Monoclonal Antibody

Catalog No	YP-Ab-04434
Isotype	IgG
Reactivity	Bovine
Applications	WB;ELISA
Gene Name	
Immunogen	Bovine serum albumin.
Specificity	BSA Monoclonal Antibody detects BSA protein.
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALB;Albumin;Bos d 6;Bos Serum albumin;Bovine Serum albumin;BSA
Observed Band	66kD
Calculated Molecular Weight	66 kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using BSA Monoclonal Antibody against BSA.