



Clenbuterol Monoclonal Antibody

Catalog No	YP-Ab-04436
Isotype	IgG
Reactivity	Species independent
Applications	ELISA
Gene Name	
Immunogen	Clenbuterol-KLH.
Specificity	Clenbuterol Monoclonal Antibody detects Clenbuterol protein.
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Clenbuterol
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
