




E7 Monoclonal Antibody

Catalog No	YP-Ab-04438
Isotype	IgG
Reactivity	Human
Applications	WB;IHC;IF;ELISA
Gene Name	
Immunogen	Purified recombinant fragment of E7 expressed in E. Coli.
Specificity	E7 Monoclonal Antibody detects E7 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	WB: 1/500 - 1/2000. IHC: 1/200 - 1/1000. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E7
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using E7 Monoclonal Antibody against recombinant protein.

Immunohistochemistry analysis of paraffin-embedded placenta tissues (left) and ovarian cancer (right) with DAB staining using E7 Monoclonal Antibody.