



HRP Monoclonal Antibody

Catalog No	YP-Ab-04442
Isotype	IgG
Reactivity	Species independent
Applications	ELISA
Gene Name	
Immunogen	Purified recombinant fragment of HRP expressed in E. Coli.
Specificity	HRP Monoclonal Antibody detects HRP protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Conjugation	HRP
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.