




Human IgG-Fc Monoclonal Antibody


Catalog No	YP-Ab-04444
Isotype	IgG
Reactivity	Human
Applications	WB;IF;ELISA
Gene Name	
Immunogen	Purified fusion protein with human IgG (Fc Specific) tag.
Specificity	Human IgG-Fc Monoclonal Antibody detects Human IgG-Fc protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COB1;YAP;YAP2;YAP65;YKI;YAP1;human IgG (Fc);Human IgG-Fc
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using Human IgG-Fc Monoclonal Antibody against different fusion proteins with human IgG (Fc specific) tag.



Confocal immunofluorescence analysis of HEK293 cells transfected with recombinant plasmid with human IgG Fc tag using Human IgG-Fc Monoclonal Antibody (green). Blue: DRAQ5 fluorescent DNA dye.