



Human IgM Monoclonal Antibody

Catalog No	YP-Ab-04446
Isotype	IgG
Reactivity	Human
Applications	ELISA
Gene Name	
Immunogen	ChromPure Human IgM (myeloma), whole molecule.
Specificity	Human IgM Monoclonal Antibody detects Human IgM protein.
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IGHM;IgM;2D10B10;MU;VH;Human IgM
Observed Band	75 kDa
Calculated Molecular Weight	69 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

