




Influenza A virus Nucleoprotein Monoclonal Antibody

Catalog No	YP-Ab-04448
Isotype	IgG
Reactivity	A/Korea/426/68(H2N2)
Applications	WB;ELISA
Gene Name	
Immunogen	Purified recombinant fragment of Influenza A virus nucleoprotein expressed in E. Coli strain BL21(DE3)
Specificity	Influenza A virus Nucleoprotein Monoclonal Antibody detects Influenza A virus Nucleoprotein protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Influenza A virus Nucleoprotein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using Influenza A virus Nucleoprotein Monoclonal Antibody against full-length recombinant Influenza A virus Nucleoprotein.