




Influenza B virus Nucleoprotein Monoclonal Antibody

Catalog No	YP-Ab-04449
Isotype	IgG
Reactivity	Stain:B/Lee/40
Applications	WB;ELISA
Gene Name	
Immunogen	Purified recombinant fragment of Influenza B virus Nucleoprotein (stain:B/Lee/40) expressed in E. Coli.
Specificity	This antibody detects endogenous levels of Influenza B virus Nucleoprotein Monoclonal Antibody at Stain:B/Lee/40, It doesn't react with total protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Influenza B virus Nucleoprotein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using Influenza B virus Nucleoprotein Monoclonal Antibody against full-length recombinant Influenza B virus Nucleoprotein.