




KSHV K8α Monoclonal Antibody


Catalog No	YP-Ab-04450
Isotype	IgG
Reactivity	KsheepV
Applications	WB;ICC;ELISA
Gene Name	
Immunogen	Purified recombinant fragment of KSHV K8α expressed in E. Coli.
Specificity	KSHV K8α Monoclonal Antibody detects KSHV K8α protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. Immunocytochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KSHV K8α
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Western Blot analysis using KSHV K8 α Monoclonal Antibody against BCBL-1 (1) and TPA induced BCBL-1 (2) cell lysate.



A



B

Immunocytochemistry analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) using KSHV K8 α Monoclonal Antibody with AEC staining.