




PAR1 Monoclonal Antibody

Catalog No	YP-Ab-04455
Isotype	IgG
Reactivity	Human
Applications	WB;ELISA
Gene Name	
Immunogen	Purified recombinant fragment of PAR1 expressed in E. Coli.
Specificity	PAR1 Monoclonal Antibody detects PAR1 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAR1;Thrombin receptor;CF2R;F2R;HTR
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using PAR1 Monoclonal Antibody against truncated GST-PAR1 recombinant protein (1).