



# Ractopamine Monoclonal Antibody

|                                    |                                                                                                                                   |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-04456                                                                                                                       |
| <b>Isotype</b>                     | IgG                                                                                                                               |
| <b>Reactivity</b>                  | Species independent                                                                                                               |
| <b>Applications</b>                | ELISA                                                                                                                             |
| <b>Gene Name</b>                   |                                                                                                                                   |
| <b>Immunogen</b>                   | Rac-OVA.                                                                                                                          |
| <b>Specificity</b>                 | Ractopamine Monoclonal Antibody detects Ractopamine protein.                                                                      |
| <b>Formulation</b>                 | Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.                                                              |
| <b>Source</b>                      | Monoclonal, Mouse                                                                                                                 |
| <b>Purification</b>                | Affinity purification                                                                                                             |
| <b>Dilution</b>                    | ELISA: 1/10000. Not yet tested in other applications.                                                                             |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                           |
| <b>Purity</b>                      | ≥90%                                                                                                                              |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                      |
| <b>Synonyms</b>                    | Ractopamine                                                                                                                       |
| <b>Observed Band</b>               |                                                                                                                                   |
| <b>Calculated Molecular Weight</b> |                                                                                                                                   |
| <b>Cell Pathway</b>                |                                                                                                                                   |
| <b>Tissue Specificity</b>          |                                                                                                                                   |
| <b>Function</b>                    |                                                                                                                                   |
| <b>Background</b>                  |                                                                                                                                   |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                              |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images

|  |
|--|
|  |
|  |
|  |
|  |
|  |
|  |
|  |