




SARS-E2 Monoclonal Antibody

Catalog No	YP-Ab-04457
Isotype	IgG
Reactivity	Human virus
Applications	WB;ELISA
Gene Name	
Immunogen	Purified recombinant fragment of SARS-E2 glycoprotein precursor expressed in E. Coli.
Specificity	SARS-E2 Monoclonal Antibody detects SARS-E2 protein.
Formulation	Subclonal supernatant.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SARS-E2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using SARS-E2 Monoclonal Antibody against SARS-E2GP3 recombinant protein.