




Human IgG mouse mAb


Catalog No	YP-Ab-04488
Isotype	IgG
Reactivity	Human
Applications	ELISA
Gene Name	human igg
Immunogen	Purified recombinant human Human IgG protein fragments expressed in E.coli.
Specificity	This antibody detects the heavy chain of Human IgG and does not cross-react with Human IgM.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	ELISA 1:10000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Human IgG;FLJ39988;FLJ40036;FLJ40253;FLJ40587;FLJ40789;FLJ40834; G1m marker;G2m marker;G3m marker;G4m marker;HDC;Heavy chain disease protein;Human immunoglobulin G;Ig gamma 1 chain C region;Ig gamma 2 chain C region;Ig gamma 3 chain C region.
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by Anti-Anthrax Protective Antigen Human antibody and then incubated with Anti-Human IgG(201068) Mouse mAb, HRP conjugated goat anti Mouse IgG was used for coloring.



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by different dilutions of supernatants which includes secreted Anti-Anthrax Protective Antigen Human antibody, HRP conjugated Anti-Human IgG(201068) Mouse mAb was used for coloring.