




HptII mouse mAb

Catalog No	YP-Ab-04490
Isotype	IgG
Reactivity	Transgenic
Applications	WB
Gene Name	hptii
Immunogen	
Specificity	Transgenic Only.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	wb dilution 1:1000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HptII
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot detection of HptII in Transgenic rice (leaf) and Non-transgenic rice (leaf) cell lysates using HptII mouse mAb (1:1000 diluted). Predicted band size: 36, 38 kDa. Observed band size: 36, 38 kDa. Beta-actin mouse mAb (200068-8F10, 1:1000 diluted) was used as a loading control.