




Phospho-Tyrosine mouse mAb

Catalog No	YP-Ab-04493
Isotype	IgG
Reactivity	Human;Mouse;Rat;AllSpeciesExpected
Applications	WB
Gene Name	py
Immunogen	Phospho-tyrosine containing peptides.
Specificity	This antibody detects endogenous levels of Phospho-Tyrosine.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	wb dilution 1:1000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Phospho-Tyrosine;py
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of extracts from A549 cells untreated and treated with Epidermal Growth Factor (10ng/ml, 15min), using Phospho-Tyrosine Mouse mAb. The phospho-specificity of the antibody was further characterized by treating the membrane with alkaline phosphatase after Western transfer.