



# GFP-Tag Polyclonal Antibody, HRP Conjugated

|                             |                                                                                                                                            |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-Ab-04572                                                                                                                                |
| Isotype                     | IgG                                                                                                                                        |
| Reactivity                  | Species independent                                                                                                                        |
| Applications                | WB;IF;                                                                                                                                     |
| Gene Name                   |                                                                                                                                            |
| Immunogen                   |                                                                                                                                            |
| Specificity                 | GFP-Tag Polyclonal Antibody HRP Conjugated, specially designed for your Western blot analysis.                                             |
| Formulation                 | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.                                                     |
| Source                      | Polyclonal, Rabbit,IgG                                                                                                                     |
| Purification                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.                              |
| Dilution                    | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000. |
| Concentration               | 1 mg/ml                                                                                                                                    |
| Purity                      | ≥90%                                                                                                                                       |
| Storage Stability           | -20°C/1 year                                                                                                                               |
| Synonyms                    | GFP-Tag;GFP Tag;GFP                                                                                                                        |
| Observed Band               |                                                                                                                                            |
| Calculated Molecular Weight |                                                                                                                                            |
| Cell Pathway                |                                                                                                                                            |
| Tissue Specificity          |                                                                                                                                            |
| Function                    |                                                                                                                                            |
| Background                  |                                                                                                                                            |
| matters needing attention   | Avoid repeated freezing and thawing!                                                                                                       |
| Usage suggestions           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.          |



## Products Images

|  |
|--|
|  |
|  |
|  |
|  |
|  |
|  |
|  |