



His-Tag Polyclonal Antibody, HRP Conjugated

Catalog No	YP-Ab-04595
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	His-Tag
Immunogen	
Specificity	His-Tag Polyclonal AntibodyFeatured Anti-His Tag Rabbit Polyclonal Antibody, HRP Conjugated, specially designed for your Western blot analysis.
Formulation	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	His-Tag;His Tag;His
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
