



Flag-Tag Monoclonal Antibody(2C5), HRP Conjugated

Catalog No	YP-Ab-04609
Isotype	IgG
Reactivity	Species independent
Applications	WB;IP;IF/ICC
Gene Name	Flag Tag
Immunogen	
Specificity	Flag-Tag Monoclonal Antibody(2C5) specially designed for your tagged protein quantification and localization.
Formulation	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Source	Monoclonal, Mouse IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000-10000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Flag Tag;Flag-Tag DDDDK TAG, DDDDK-TAG;Flag
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
