



Myc-Tag Monoclonal Antibody(3E8), AbFluor™ 594 Conjugated

Catalog No	YP-Ab-04635
Isotype	IgG
Reactivity	Species independent
Applications	WB;ELISA;IP;IF
Gene Name	Myc Tag
Immunogen	
Specificity	Myc-Tag Monoclonal Antibody(3E8) AbFluor™ 594 Conjugated specially designed for your Immunofluorescence analysis.
Formulation	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Source	Monoclonal, Mouse IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF: 1:1000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Myc-Tag;Myc Tag;Myc
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images