



# T7-Tag Monoclonal Antibody(1H5), HRP Conjugated

|                                    |                                                                                                                                            |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-04662                                                                                                                                |
| <b>Isotype</b>                     | IgG                                                                                                                                        |
| <b>Reactivity</b>                  | Species independent                                                                                                                        |
| <b>Applications</b>                | WB                                                                                                                                         |
| <b>Gene Name</b>                   |                                                                                                                                            |
| <b>Immunogen</b>                   |                                                                                                                                            |
| <b>Specificity</b>                 | T7-Tag Monoclonal Antibody(1H5) HRP Conjugated, specially designed for your Western blot analysis.                                         |
| <b>Formulation</b>                 | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.                                                     |
| <b>Source</b>                      | Monoclonal, Mouse IgG                                                                                                                      |
| <b>Purification</b>                | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.                                 |
| <b>Dilution</b>                    | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000. |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                    |
| <b>Purity</b>                      | ≥90%                                                                                                                                       |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                               |
| <b>Synonyms</b>                    | T7-Tag;T7                                                                                                                                  |
| <b>Observed Band</b>               |                                                                                                                                            |
| <b>Calculated Molecular Weight</b> |                                                                                                                                            |
| <b>Cell Pathway</b>                |                                                                                                                                            |
| <b>Tissue Specificity</b>          |                                                                                                                                            |
| <b>Function</b>                    |                                                                                                                                            |
| <b>Background</b>                  |                                                                                                                                            |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                       |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.          |



## Products Images

|  |
|--|
|  |
|  |
|  |
|  |
|  |
|  |
|  |