



# V5-Tag Monoclonal Antibody(3C8), AbFluor™ 488 Conjugated

|                             |                                                                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-Ab-04668                                                                                                                                 |
| Isotype                     | IgG                                                                                                                                         |
| Reactivity                  | Species independent                                                                                                                         |
| Applications                | WB;IP;IF                                                                                                                                    |
| Gene Name                   |                                                                                                                                             |
| Immunogen                   |                                                                                                                                             |
| Specificity                 | V5-Tag Monoclonal Antibody(3C8) AbFluor™ 488 Conjugated specially designed for your Immunofluorescence analysis.                            |
| Formulation                 | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.                                                      |
| Source                      | Monoclonal, Mouse IgG                                                                                                                       |
| Purification                | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.                                  |
| Dilution                    | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF: 1:1000. |
| Concentration               | 1 mg/ml                                                                                                                                     |
| Purity                      | ≥90%                                                                                                                                        |
| Storage Stability           | -20°C/1 year                                                                                                                                |
| Synonyms                    | V5-Tag;V5                                                                                                                                   |
| Observed Band               |                                                                                                                                             |
| Calculated Molecular Weight |                                                                                                                                             |
| Cell Pathway                |                                                                                                                                             |
| Tissue Specificity          |                                                                                                                                             |
| Function                    |                                                                                                                                             |
| Background                  |                                                                                                                                             |
| matters needing attention   | Avoid repeated freezing and thawing!                                                                                                        |

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**