



VSV-G-Tag Monoclonal Antibody(8D6), AbFluor™ 680 Conjugated

Catalog No	YP-Ab-04683
Isotype	IgG
Reactivity	Species independent
Applications	WB;IP;IF
Gene Name	
Immunogen	A synthetic peptide corresponding to aa 501-511 of vesicular stomatitis virus glycoprotein (VSV-G). This vector is a pseudo lentiviral vector.
Specificity	VSV-G-Tag Monoclonal Antibody(8D6) AbFluor™ 680 Conjugated specially designed for your Immunofluorescence analysis.
Formulation	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Source	Monoclonal, Mouse IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF: 1:1000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VSV-G-Tag
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images