




Strep-Tag Monoclonal Antibody(7A10)

Catalog No	YP-Ab-04738
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	
Immunogen	Synthetic Peptide of Strep-Tag
Specificity	The antibody detects Strep tag fusion proteins.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB: 1:5000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Strep;Strep-Tag;Strep-Tag II
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



1ug Strep II fusion protein+ Primary antibody dilution at
1) 1:5000 2) 1:10000