



GFP Monoclonal Antibody(Mix)

Catalog No	YP-Ab-04749
Isotype	IgG
Reactivity	Species independent
Applications	WB;IF/ICC;IP
Gene Name	
Immunogen	Recombinant Protein of GFP
Specificity	The antibody detectsThe antibody detects RFP, YFP, CFP.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB: 1:5000 IP: 1:200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	27KDa
Calculated Molecular Weight	27KD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
