



mStrawberry Monoclonal Antibody(Mix)

Catalog No	YP-Ab-04759
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	
Immunogen	Synthetic Peptide of mStrawberry
Specificity	The antibody detects mStrawberry & mStrawberry Tag protein. It May had cross reactivity with mRFP protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB: 1:2000-5000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	mStrawberry
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	Living Colors mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP. mStrawberry excitation and emission maxima are 574 and 596 nm, respectively.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of mStrawberry recombinant protein, diluted at 1:5000.