



mScarlet Mouse mAb(Mix-mA)


Catalog No	YP-Ab-04830
Isotype	IgG
Reactivity	All
Applications	IHC;WB
Gene Name	
Immunogen	Synthesized peptide derived from mScarlet
Specificity	This antibody detects endogenous levels of mScarlet
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.72% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p1:50-200 ,WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	mScarlet(Mix-mA);mScarlet
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	mScarlet: a bright monomeric red fluorescent protein for cellular imaging.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

94KD —————
66KD —————
45KD —————
35KD —————
27KD —————

14KD —————



Western blot analysis of mScarlet recombinant transfected 293T Cell Lysate using mScarlet Mouse mAb diluted at 1:5,000.