



Digoxin Monoclonal Antibody(4F6)

Catalog No	YP-Ab-04841
Isotype	IgG
Reactivity	0
Applications	ELISA
Gene Name	
Immunogen	Digoxin
Specificity	
Formulation	
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	ELISA 1/10000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	283-75-5; CasNo 283-75-5;  ; Digoxin
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	The most common indications for digoxin are atrial fibrillation and atrial flutter with rapid ventricular response, though beta blockers and/or calcium channel blockers are a better first choice. There is tentative evidence that digoxin may increase the risk of death, though another meta-analysis reports no change in mortality. High ventricular rate leads to insufficient diastolic filling time. By slowing down the conduction in the AV node and increasing its refractory period, digoxin can reduce the ventricular rate. The arrhythmia itself is not affected, but the pumping function of the heart improves, owing to improved filling.

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images