



Glu-Glu Tag mouse Monoclonal Antibody(2G1)

Catalog No	YP-Ab-04874
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	
Immunogen	Synthetic Peptide of Glu-Glu Tag
Specificity	The antibody detects Glu-Glu Tag protein
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	YPYDVPDYA;Glu-Glu Tag
Observed Band	
Calculated Molecular Weight	1 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Epitope tags are useful for the labeling and detection of proteins using immunoblotting, immunoprecipitation, and immunostaining techniques. Because of their small size, they are unlikely to affect the tagged protein's biochemical properties.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images