




Placental Alkaline Phosphatase(PLAP) mouse mAb(PT2105)

Catalog No	YP-Ab-04886
Isotype	IgG
Reactivity	Human
Applications	IHC;IF;WB
Gene Name	PLAP
Immunogen	Synthesized peptide derived from human Placental Alkaline Phosphatase(PLAP)
Specificity	This antibody detects endogenous levels of human Placental Alkaline Phosphatase(PLAP)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-300, WB 1:500-2000. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PLAP-1;PALP;EC:3.1.3.1;ALPP;ALP;Placental Alkaline Phosphatase(PLAP) (PT2105);PLAP;Placental Alkaline Phosphatase(PLAP)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human Seminoma Antibody was diluted at 1:200(4° overnight).