




AHNK Polyclonal Antibody

Catalog No	YP-Ab-04887
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	AHNAK PM227
Immunogen	Synthesized peptide derived from human protein . at AA range: 5720-5800
Specificity	AHNK Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC-p 1:50-300. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PM227;Neuroblast differentiation-associated protein AHNAK;Desmoyokin;AHNK;AHNAK PM227;Neuroblast differentiation-associated protein AHNAK (Desmoyokin)
Observed Band	647kD
Calculated Molecular Weight	647kD
Cell Pathway	Nucleus.
Tissue Specificity	Cervix,Epithelium,Kidney,Liver,PCR rescued clones,Placenta,Tongue,Uterus,
Function	function:May be required for neuronal cell differentiation.,similarity:Contains 1 PDZ (DHR) domain.,subunit:Interacts with DYSF; the interaction is direct and Ca(2+)-independent.,
Background	function:May be required for neuronal cell differentiation.,similarity:Contains 1 PDZ (DHR) domain.,subunit:Interacts with DYSF; the interaction is direct and Ca(2+)-independent.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Tris-EDTA, pH9.0 was used for antigen retrieval. 2 Antibody was diluted at 1:200 (4° overnight). 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).