



MSX1 Polyclonal Antibody

Catalog No	YP-Ab-04969
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	MSX1 HOX7
Immunogen	Synthesized peptide derived from human protein . at AA range: 70-150
Specificity	MSX1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Homeobox protein Hox 7;Homeobox protein MSX 1;HOX7;HYD1;msh homeobox 1;Msh;MSX1;MSX1 HOX7;Homeobox protein MSX-1 (Homeobox protein Hox-7) (Msh homeobox 1-like protein)
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	Nucleus.
Tissue Specificity	Expressed in the developing nail bed mesenchyme.
Function	disease:A chromosomal aberration involving MSX1 is a cause of Wolf-Hirschhorn syndrome (WHS) [MIM:194190]. WHS is caused by sub-telomeric deletions in the short arm of chromosome 4. WHS is characterized by profound mental retardation, heart defects, and facial clefting.,disease:Defects in MSX1 are a cause of autosomal dominant hypodontia (HYD1) [MIM:106600]; also known as familial or selective tooth agenesis. Absence of less than 6 teeth is referred to as hypodontia. Agenesis of one or more teeth constitutes one of the most common developmental anomalies in man. Reported incidences vary from 1.6% to 9.6%, excluding third molar (Wisdom tooth) agenesis, which occurs in 20% of the population.,disease:Defects in MSX1 are the cause of non-syndromic orofacial cleft type 5 (OFC5) [MIM:608874]; also called non-



syndromic cleft lip with or without cleft palate 5. Non-syndromic orofacial cleft is a

Background

This gene encodes a member of the muscle segment homeobox gene family. The encoded protein functions as a transcriptional repressor during embryogenesis through interactions with components of the core transcription complex and other homeoproteins. It may also have roles in limb-pattern formation, craniofacial development, particularly odontogenesis, and tumor growth inhibition. Mutations in this gene, which was once known as homeobox 7, have been associated with nonsyndromic cleft lip with or without cleft palate 5, Witkop syndrome, Wolf-Hirschorn syndrome, and autosomal dominant hypodontia. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night