



DNJB5 Polyclonal Antibody


Catalog No	YP-Ab-05068
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	DNAJB5 HSC40
Immunogen	Synthesized peptide derived from human protein . at AA range: 140-220
Specificity	DNJB5 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DNJB5;DNAJB5 HSC40;DnaJ homolog subfamily B member 5 (Heat shock protein Hsp40-2) (Heat shock protein Hsp40-3) (Heat shock protein cognate 40) (Hsc40)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	nucleus,cytosol,
Tissue Specificity	Brain,Colon,Skin,
Function	induction:Expressed under normal conditions, its expression can further be increased after various stress treatments.,similarity:Contains 1 J domain.,
Background	DNAJB5 belongs to the evolutionarily conserved DNAJ/HSP40 protein family. For background information on the DNAJ family, see MIM 608375.[supplied by OMIM, Mar 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC12 cells,
primary antibody was diluted at 1:1000, 4° over night