



IFI44 Polyclonal Antibody

Catalog No	YP-Ab-05085
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	IFI44 MTAP44
Immunogen	Synthesized peptide derived from human protein . at AA range: 310-390
Specificity	IFI44 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IFI44;interferon induced protein 44;Interferon-induced protein 44;Microtubule-associated protein 44;MTAP44;IFI44 MTAP44;Interferon-induced protein 44 (p44) (Microtubule-associated protein 44)
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Lung,
Function	function:This protein aggregates to form microtubular structures.,induction:By alpha and beta interferons, but not by gamma interferons.,similarity:Belongs to the IFI44 family.,
Background	function:This protein aggregates to form microtubular structures.,induction:By alpha and beta interferons, but not by gamma interferons.,similarity:Belongs to the IFI44 family.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images