



IFIT2 Polyclonal Antibody

Catalog No	YP-Ab-05086
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	IFIT2 CIG-42 G10P2 IFI54 ISG54
Immunogen	Synthesized peptide derived from human protein . at AA range: 200-280
Specificity	IFIT2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	cig42;CIG-42;G10P2;GARG 39;IFI 54;IFIT2;IFIT2 CIG-42 G10P2 IFI54 ISG54; Interferon-induced protein with tetratricopeptide repeats 2 (IFIT-2) (ISG-54 K) (Interferon-induced 54 kDa protein) (IFI-54K) (P54)
Observed Band	50-55 kDa
Calculated Molecular Weight	484 aa, 56 kDa
Cell Pathway	Cytoplasm . Endoplasmic reticulum .
Tissue Specificity	Bone marrow,Epithelium,Testis,
Function	induction:By interferons.,similarity:Belongs to the IFIT family.,similarity:Contains 6 TPR repeats.,
Background	induction:By interferons.,similarity:Belongs to the IFIT family.,similarity:Contains 6 TPR repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
