



IFIT3 Polyclonal Antibody

Catalog No	YP-Ab-05087
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	IFIT3 CIG-49 IFI60 IFIT4 ISG60
Immunogen	Synthesized peptide derived from human protein . at AA range: 330-410
Specificity	IFIT3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CIG 49;CIG49;CIG-49;GARG 49;GARG49;IFIT3;IFIT3 CIG-49 IFI60 IFIT4 ISG60;Interferon-induced protein with tetratricopeptide repeats 3 (IFIT-3) (CIG49) (ISG-60) (Interferon-induced 60 kDa protein) (IFI-60K) (Interferon-induced protein with tetratricopeptide repeats 4) (IFIT-
Observed Band	60 kDa
Calculated Molecular Weight	56 kDa
Cell Pathway	Cytoplasm . Mitochondrion .
Tissue Specificity	Expression significantly higher in peripheral blood mononuclear cells (PBMCs) and monocytes from systemic lupus erythematosus (SLE) patients than in those from healthy individuals (at protein level). Spleen, lung, leukocytes, lymph nodes, placenta, bone marrow and fetal liver.
Function	similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,
Background	similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images