




IFT1B Polyclonal Antibody

Catalog No	YP-Ab-05100
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	IFIT1B IFIT1L
Immunogen	Synthesized peptide derived from human protein . at AA range: 10-90
Specificity	IFT1B Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IFT1B;IFIT1B IFIT1L;Interferon-induced protein with tetratricopeptide repeats 1B (Interferon-induced protein with tetratricopeptide repeats 1-like protein)
Observed Band	52kD
Calculated Molecular Weight	52kD
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,
Background	similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000