



IIGP5 Polyclonal Antibody

Catalog No	YP-Ab-05102
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	IRGC IIGP5 IRGC1
Immunogen	Synthesized peptide derived from human protein . at AA range: 350-430
Specificity	IIGP5 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IRGC1;Interferon-inducible GTPase 5;Immunity-related GTPase cinema 1;IIGP5; EC:3.6.5.-;IRGC IIGP5 IRGC1;Interferon-inducible GTPase 5 (EC 3.6.5.-) (Immunity-related GTPase cinema 1)
Observed Band	50 kDa
Calculated Molecular Weight	50 kDa
Cell Pathway	membrane,
Tissue Specificity	Brain,Testis,
Function	catalytic activity:GTP + H(2)O = GDP + phosphate.,similarity:Belongs to the interferon-inducible GTPase family.,
Background	catalytic activity:GTP + H(2)O = GDP + phosphate.,similarity:Belongs to the interferon-inducible GTPase family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
