



ITFG2 Polyclonal Antibody

Catalog No	YP-Ab-05112
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ITFG2
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	ITFG2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Integrin-alpha FG-GAP repeat-containing protein 2;KICSTOR complex protein ITFG2;MDS028;ITFG2
Observed Band	52 kDa
Calculated Molecular Weight	49 kDa
Cell Pathway	Lysosome membrane . Localization to lysosomes is amino acid-independent. .
Tissue Specificity	Hematopoietic stem/progenitor cells,Ovary,
Function	similarity:Contains 2 FG-GAP repeats.,
Background	similarity:Contains 2 FG-GAP repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
