




MAGF1 Polyclonal Antibody

Catalog No	YP-Ab-05135
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	MAGEF1
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	MAGF1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MAGF1;MAGEF1;Melanoma-associated antigen F1 (MAGE-F1 antigen)
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	extracellular exosome,
Tissue Specificity	Ubiquitous.
Function	similarity:Contains 1 MAGE domain.,tissue specificity:Ubiquitous.,
Background	This intronless gene encodes a member of the MAGE superfamily. It is ubiquitously expressed in normal tissues and in tumor cells. This gene includes a microsatellite repeat in the coding region. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4°over night