



PDCD7 Polyclonal Antibody

Catalog No	YP-Ab-05148
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	PDCD7
Immunogen	Synthesized peptide derived from human protein . at AA range: 300-380
Specificity	PDCD7 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ES18;HES18;PDCD7;programmed cell death 7;Programmed cell death protein 7 (ES18) (hES18)
Observed Band	55 kDa
Calculated Molecular Weight	485 aa, 55 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Duodenum,Placenta,Small intestine,Thymus,
Function	function:Promotes apoptosis when overexpressed.,subunit:Interacts with RBM40.,
Background	This gene encodes a 59 kDa protein that is associated with the U11 small nuclear ribonucleoprotein (snRNP), which is a component of the minor U12-type spliceosome responsible for catalyzing pre-mRNA splicing of U12-type introns. [provided by RefSeq, Dec 2010],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images