




PDD2L Polyclonal Antibody

Catalog No	YP-Ab-05149
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	PDCD2L
Immunogen	Synthesized peptide derived from human protein . at AA range: 180-260
Specificity	PDD2L Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PDD2L;PDCD2L;Programmed cell death protein 2-like
Observed Band	39kD
Calculated Molecular Weight	39kD
Cell Pathway	cytoplasm,membrane,
Tissue Specificity	Higher expression in lung, colon, mammary gland, cervix, stomach and small intestine.
Function	
Background	PDCD2L (Programmed Cell Death 2 Like) is a Protein Coding gene. An important paralog of this gene is ENSG00000266953. Over-expression suppresses AP1, CREB, NFAT, and NF-κB transcriptional activation, and delays cell cycle progression at S phase.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4°over night