



STRN4 Polyclonal Antibody

Catalog No	YP-Ab-05188
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	STRN4 ZIN
Immunogen	Synthesized peptide derived from human protein . at AA range: 250-330
Specificity	STRN4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Striatin-4;ZIN;Zinedin;STRN4;STRN4 ZIN;Striatin-4 (Zinedin)
Observed Band	82kD
Calculated Molecular Weight	82kD
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein . Cell projection, dendritic spine . CTTNBP2-binding may regulate dendritic spine distribution. .
Tissue Specificity	Brain,Muscle,Testis,
Function	function:Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).,similarity:Belongs to the WD repeat striatin family.,similarity:Contains 7 WD repeats.,
Background	function:Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).,similarity:Belongs to the WD repeat striatin family.,similarity:Contains 7 WD repeats.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000