



# BCL6B Polyclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-05198                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Reactivity</b>                  | Human;Mouse                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Applications</b>                | WB;ELISA                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Gene Name</b>                   | BCL6B BAZF ZNF62                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen</b>                   | Synthesized peptide derived from human protein . at AA range: 100-180                                                                                                                                                                                                                                                                                                                                                           |
| <b>Specificity</b>                 | BCL6B Polyclonal Antibody detects endogenous levels of protein.                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                                                           |
| <b>Dilution</b>                    | WB 1:500-2000 ELISA 1:5000-20000                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Synonyms</b>                    | BCL6B;BCL6B BAZF ZNF62;B-cell CLL/lymphoma 6 member B protein (Bcl6-associated zinc finger protein) (Zinc finger protein 62)                                                                                                                                                                                                                                                                                                    |
| <b>Observed Band</b>               | 52kD                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Calculated Molecular Weight</b> | 52kD                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Cell Pathway</b>                | Nucleus .                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Tissue Specificity</b>          | Ubiquitously expressed with higher expression found in heart and placenta.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Function</b>                    | function:Acts as a sequence-specific transcriptional repressor in association with BCL6. May function in a narrow stage or be related to some events in the early B-cell development.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 C2H2-type zinc fingers.,subunit:Associates with BCL6 through the BTB domain.,tissue specificity:Ubiquitously expressed with higher expression found in heart and placenta., |
| <b>Background</b>                  | function:Acts as a sequence-specific transcriptional repressor in association with BCL6. May function in a narrow stage or be related to some events in the early B-cell development.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 C2H2-type zinc fingers.,subunit:Associates with BCL6 through the BTB domain.,tissue specificity:Ubiquitously expressed with higher expression found                         |



in heart and placenta.,

**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images