



TNIP3 Polyclonal Antibody

Catalog No	YP-Ab-05207
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TNIP3 ABIN3 LIND
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	TNIP3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TNIP3;TNIP3 ABIN3 LIND;TNFAIP3-interacting protein 3 (A20-binding inhibitor of NF-kappa-B activation 3) (ABIN-3) (Listeria-induced gene protein)
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	
Tissue Specificity	Highly expressed in lung, lymph node, thymus and fetal liver. Expressed at lower levels in bone marrow, brain, kidney, spleen, leukocytes and tonsils. Could be detected in heart, salivary gland, adrenal gland, pancreas, ovary and fetal brain. High levels detected in liver, colon, small intestine, muscle, stomach, testis, placenta, thyroid, uterus, prostate, skin and PBL.
Function	induction:By Listeria infection. Expression is slightly down-regulated by dexamethasone and slightly up-regulated by IL-10.,tissue specificity:Highly expressed in lung, lymph node, thymus and fetal liver. Expressed at lower levels in brain, spleen, leukocytes and tonsils.,
Background	induction:By Listeria infection. Expression is slightly down-regulated by dexamethasone and slightly up-regulated by IL-10.,tissue specificity:Highly expressed in lung, lymph node, thymus and fetal liver. Expressed at lower levels



in brain, spleen, leukocytes and tonsils.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night