



PO3F3 Polyclonal Antibody

Catalog No	YP-Ab-05212
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	POU3F3 BRN1 OTF8
Immunogen	Synthesized peptide derived from human protein . at AA range: 90-170
Specificity	PO3F3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PO3F3;POU3F3 BRN1 OTF8;POU domain;class 3;transcription factor 3 (Brain-specific homeobox/POU domain protein 1) (Brain-1) (Brn-1) (Octamer-binding protein 8) (Oct-8) (Octamer-binding transcription factor 8) (OTF-8)
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain.
Function	similarity:Belongs to the POU transcription factor family. Class-3 subfamily., similarity:Contains 1 homeobox DNA-binding domain.,similarity:Contains 1 POU-specific domain.,tissue specificity:Brain.,
Background	This gene encodes a POU-domain containing protein that functions as a transcription factor. The encoded protein recognizes an octamer sequence in the DNA of target genes. This protein may play a role in development of the nervous system. [provided by RefSeq, Apr 2015],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images