




ABHD1 Polyclonal Antibody

Catalog No	YP-Ab-05261
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	ABHD1 LABH1
Immunogen	Synthesized peptide derived from human protein . at AA range: 220-300
Specificity	ABHD1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABHD 1;Abhydrolase domain-containing protein 1;Alpha/beta hydrolase domain-containing protein 1;EC:3.1.1.-;LABH1;ABHD1;ABHD1 LABH1;Abhydrolase domain-containing protein 1 (EC 3.1.1.-) (Lung alpha/beta hydrolase 1)
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Belongs to the AB hydrolase superfamily. AB hydrolase 4 family.,tissue specificity:Ubiquitously expressed.,
Background	This gene is a member of the AB hydrolase superfamily and encodes a protein with an alpha/beta hydrolase fold. This domain is common to a number of hydrolytic enzymes of widely differing phylogenetic origins and catalytic functions. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night